


AF-FL035 VITA EMOTIONAL TUNNEL MIX

2 Showerheads 35x15 cm and 1 showerhead 35x35 cm with different functions controlled by program station P600 for sequential control of each step. Water jets, flow rate ,chromotherapy,times including electronic mixer for control temperature. Possibility to combine liquid or dry aromatherapy and music.

VITA EMOTIONAL TUNNEL MIX




SHOWERHEAD
N3515C15
R35C12
N3515C3A

CASE
I35-3



P600



SEQUENTIAL PROGRAM STATION PLUS FOR SHOWER PATH.

Electronic power station prepared to activate in sequence 8 steps with temporary control of water for each step. Light 1 channel RGB. Complete of Electronic mixer for control temperature. Equipped with keyboard with button to start program with transformer 12 V. Optional aromatherapy during the entire program.

N.B. To be teamed with a single wiring electrovalve C320 12V for each step.

C320



VALVES.

Solenoid valve for program station P600


P02214



CABLES.

6-pole cable with M12 female quick lock connector L = 5 m to be teamed with shower system and program station

INSTALLATION VITA EMOTIONAL TUNNEL MIX



LEGEND

- WATER SUPPLY 1/2"
- ELECTRICAL CABLES
- - - - - CORRUGATED PIPE DIAM. 20mm
- - - - - LED CABLES DIAM. 20mm

TECHINCAL INFORMATION

WATER FLOW	min 15 l/m
WATER SUPPLY	1/2" inch
OPERATING PRESSURE	1,5/5 bar
RECOMMENDED PRESSURE	3 bar

CONFIGURATION

Code	Description	Jets	Quantity
EMOTIONAL MIX	SHOWERHEAD		1
P600	PROGRAM STATION	-	1
C320	VALVES	-	5
G0279	CONTROL BOX	-	1
P02214	CABLES	-	1

OPTIONAL

OPTIONAL 01 - INVERTED PATH

Code	Description	Quantity
C45I	KEYBOARD WITH OPTIONAL BUTTON Built-in keyboard in stainless steel, flush with the wall, with 1 IP67 button and 3 modules flush mounting box (type 503), to activate opposite path for sequential power station our cod. P300	1

OPTIONAL 01 - MUSIC

Code	Description	Quantity
PSOUND	mp3 device including amplifier with relative interface for sequential program station	1
KSOUND	Kit audio with 2 high output exit and connectors 2 poles . N.B. -To be teamed with keyboard touch (our cod. C57- C60) N.B. -To be teamed with PSOUND for use with aquaform program station.	6

OPTIONAL 02 - LIQUID AROMA

Code	Description	Quantity
PE427	Liquid essences kit. Dosing pump complete of: - Injector connector with tube kit - Neck off tank - Tank wall bracket	1
EL1/EL2/EL3/EL4	5 lt mint essence tank	1

OPTIONAL 03 - DRY AROMA

Code	Description	Quantity
A40	Recessed diffuser for aromatherapy system for diffusion and relative control flush-mounted 4 module box (type 504) kit packaging 3 flower-relax-energy essences supplied. Suitable for humid environments.	1
E1/E2/E3/E4/E5	AROMATHERAPY ESSENCE 3 cartridges set.	1